SILICONISED LATEX FOLEY CATHETER

INSTRUCTION FOR USE

DESCRIPTION

All catheters basically comprise a 2-lumen with proximal funnel, inflation valve and distal retaining balloon.

Balloon fill volumes are in ml. And the shaft size is in French (Fr.), Charrière (Ch.) which is indicated on the funnel of each individual catheter and on individual pack label and outer carton.

Sizes range from Ch./Fr. 6 to 26.

Length: 30 cm for sizes 6 – 10 Ch.
22 cm for female version, sizes 12 – 26 Ch.
41 cm for male and irrigation version, sizes 12-26 Ch.

The Securdrain balloon catheter is made of siliconized latex.

Dimension (Ch.)	Balloon Capacity
6	3 ml
8	3/5 ml.
10	3/5 ml
12-26	5-15 ml
12-26	30 ml

INDICATIONS

These catheters are indicated for routine drainage of the bladder or for routine post-operative drainage and irrigation bladder.

WARNING!!

THE CH 6 FOLEY'S DIMENSIONS IMPOSE, TO THE PERSONNEL IN CHARGED OF ITS MANIPULATION ON THE END-USER., TO USE CAUTION AND GENTLENESS. IN PARTICULAR THE BALLOON INFLATION/DEFLATION OPERATION MUST BE EFFECTED WITH PLENTY OF SENSITIVENESS BY EXERTING SOFT PRESSURE ON THE SYRINGE AND BY INSERTING GRADUALLY THE LIQUID.

Read all instructions carefully before use:

PREPARATION FOR USE:

Remove catheters from its sterile pouch using aseptic technique. Visually inspect the catheter for any mechanical damage. If a trial inflation of the balloon is carried out, ensure the balloon is empty prior to insertion. Only sterile water should be used to fill the balloon.

DIRECTIONS FOR USE

Lubricate the catheter using only water based lubricant or gel (with or without local anaesthetic as required). Do Not use petroleum based lubricants on catheter. Using sterile technique as indicated in published literature; insert the catheter fully into the urethra ensuring that the balloon is beyond the bladder neck.

A luer syringe may be used to inflate the balloon.

Inflate the balloon with sterile water.

The capacity of the balloon is indicated on the funnel of the catheter or on the label of the packaging.

DO NOT OVER INFLATE.

Ensure function by checking for urine flow from the drainage funnel. To deflate the balloon prior to removal, insert a Luer syringe as appropriate. Aspirate the syringe gently to remove the inflation fluid. DO not use excessive aspiration on the syringe during deflation as this may cause a vacuum collapse of the inflation lumen, which may impair aspiration of inflation fluid.

Maximum recommended indwelling time is 7 (seven) days.

The inflation balloon holds the catheters in the bladder. If the balloon loses fluid, e.g. by diffusion fluid through the balloon membrane, the catheter can, with time, slide out the bladder accidentally.

The patient should be routinely monitored in accordance with accepted procedures and the catheter removed after a suitable interval, no longer than 7 days, as determined by a physician or other suitably qualified personnel.


POTENTIAL COMPLICATIONS

If difficulty is encountered aspirating the balloon with a syringe, a rare and infrequent reported event, it is possible to use one of the following procedures:


- 1. Deflation through torsion:** aspire drawing slowly the catheter extremity, by effecting rotation movement alternatively towards left and right in order to remove any possible obstructing entity from the drainage channel of the balloon.


- 2. Cut of the distal extremity of the catheter:** because of the frequent torsions, on the final part of the catheter there are often encrustations. In this case, it is possible - once verified the site of the stenosis and/or encrustation - to cut the distal part of the catheter. For security reason, it is advisable to maintain a distance of about 5 cm from the urethral meatus.


- 3. Introduction of a tube and perforation:** introduce a little tube in the balloon inflation channel (for example, the urethral catheter mandrel of the mandrel of a catheters for angiography). Try then to perforate the balloon with the inserted mandrel (**warning!** The balloon do not have to burst! Then verify the complete removal).


- 4. Percutaneous punching:** this method can be used both through a radiological exam and through an ultrasound exam. During the punching operation, the balloon has to be positioned inside the vesical neck through a gentle traction on the catheter. As punching tool can be utilised needles for nephrostomy, for biopsy or indwelling catheter cannulae (**warning!** The balloon do not have to burst! Then verify the complete removal).

- 5. Harpoon technique:** after having cut the catheter extremity, introduce over the catheter a cytoscope previously coated with lubricating gel. Lead the cytoscope to the balloon. Then insert inside the cytoscope a mandrel in order to punch the balloon using the O'Flynn technique (**warning!** The balloon do not have to burst! Then verify the complete removal).

- 6. transrectal or transvaginal punching:** by using a flexible biopsy needle, the balloon can be perforated punching inside the vesical neck (**warning!** The balloon do not have to burst! Then verify the complete removal).

Incorrectly positioned catheters can cause urethral damage if the balloon is inflated within the urethra. Irritation of the urethra mucosa, blockage of the catheters due to encrustation and catheter-induced infections are documented complications with some patients.

CAUTION

For urological use only. Use only with water based lubricant. Always inflate the balloon with sterile water.

Never clamp the catheter: a catheter plug or valve may be used if required.

Never use a syringe to puncture the catheter shaft for urine sampling. There is a danger of needle stick injury and the catheter may be compromised. Phthalates (DEHP) free.

PRESENTATION

Each product is supplied in a individual sterile pouch. The product is guaranteed sterile unless the individual packaging has been opened or damaged.

This product is for single use only. Risk of infections in case of reuse.

DO NOT RESTERILIZE OR REUSE.

Dispose of according to applicable regulation.

- CATETERE TIPO FOLEY IN LATTICE SILICONATO
- SILICONIZED LATEX FOLEY TYPE CATHETER

STERILE EO

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Review date: 21/10/2014

and Importer:
medical supplies and accessories
Securmed
SECURMED S.p.A.
VIA MONTE GRAPPA, 2/G
36016 THIENE (VI) - ITALY

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