

NOBAMED Paul Danz AG

# **NOBAHAFT®-fein**

### **REF 470206**

## Product Description, Intended use, Application

The elastic fixation bandage NOBAHAFT<sup>®</sup>fein, white, 6 cm x 20 m, is latex-free and cohesive. It does not adhere to the skin or clothing, but only to itself. The bandage has a fine, smooth structure. It is used under tape bandages, mainly at the ankle joint, but may also be used to fix wound pads.

### Composition

Viscose, polyamide, latex-free polymer

### Contraindications

The product should not be used in the case of a known allergy against the material.

### Note

During application, cohesive bandages should only be stretched slightly as otherwise the arterial blood supply might be disturbed.

#### **Incident reporting**

According to MDR (EU) 2017/745, if serious incidents occur in relation to the device, they must be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

Normative and Regulative Requirements, Common Standards Medical Device according to MDD 93/42/EEC, MDR (EU) 2017/745. The product does not contain dangerous toxic substances according to REACH.

### Packaging

Primary packaging:

folding box made of cellulose

Secondary packaging: carton made of cellulose

Symbols used in labelling Explanation at www.nobamed.com



Marking on all packaging levels with CE and according to DIN EN ISO 15223-1and DIN EN ISO 20417

**Storage and Transport** Dry and dustfree

#### Single use device

Reusing a single use medical device can lead to microbiological danger. Reprocessing for reuse can decrease the product's performance significantly.

#### Disposal

According to locally applicable legal regulations and standards of infection prophylaxis.